

Optimal Twillkudde

Av högkvalitativ bollfiber

Optimal Twillkudde har en fyllning av bollfiber i syfte att uppnå en god komfort.

Optimal Twillkudde är utvecklad för vårdens behov och passar även i hemmiljö.

Överdrag för god hygien och flexibilitet

Som tillbehör finns ett vätsketätt hygienöverdrag. Det vätsketäta hygienöverdraget är tillverkat i ett tyg med fyrvägsstretch⁽¹⁾⁽²⁾ för att minska risken för skjuvkrafter, samt är ånggenomsläppligt⁽³⁾ för att reducera risken för hudmaceration. Hygienöverdraget har även en vätsketät dragkedja.



Teknisk specifikation	
Fyllning	Polyester/bollfiber
Storlekar	30x45 cm, 50x60 cm, 50x70 cm
Tillgängliga överdrag	Hygienöverdrag (polyuretanbelagd polyester)
Instruktion för rengöring	Rengöring av överdrag: Avtorkning med rengörings- och/eller desinfektionsmedel. Maskintvätt twillkudde max 70 °C
CE-märkning	Registrerad och CE-märkt i enlighet med MDR (EU) 2017/745

Referenser

- (1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.
- (2) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements.
- (3) SS-EN ISO 15496:2004, DIN 53122-1.



Läs alltid bruksanvisningen innan användning.

Care of Sweden and Curera are trademarks of Care of Sweden AB. © Care of Sweden AB, 2024. All rights reserved. Last update: 2024-01-22.