

Optimal Hygientopp

Slitstarkt madrasskydd med bra passform & följsamhet

Optimal Hygientopp är ett enkelt, slitstarkt och funktionellt madrasskydd med bra passform och följsamhet. Optimal Hygientopp är ett avtagbart och vätsketätt hygienöverdrag för enkel rengöring.

Det vätsketäta hygienöverdraget är tillverkat i ett tyg med fyrvägsstretch⁽¹⁾⁽²⁾ för att minska risken för skjuvkrafter, samt är ånggenomsläppligt⁽³⁾ för att reducera risken för hudmaceration.



Teknisk specifikation	
Material	Polyuretan/polyester
Storlekar	80x205 cm 85x200 cm 90x200 cm 105x200 cm 120x200 cm
Instruktion för rengöring	Avtorkning med rengöringsoch/ eller desinfektionsmedel. Maskintvätt max 95 °C, torktumling.
CE-märkning	Registrerad och CE-märkt i enlighet med MDR (EU) 2017/745

Referenser

- (1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.
- (2) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements.
- (3) SS-EN ISO 15496:2004, DIN 53122-1.



Läs alltid bruksanvisningen innan användning.

Care of Sweden and Curera are trademarks of Care of Sweden AB. © Care of Sweden AB, 2024. All rights reserved. Last update: 2024-01-22.